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Acc. to ISO 8637-1; *in vitro* data are likely to differ from *in vivo* data due to the patient's blood composition and clinical settings.

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Blood connectors are compatible with ISO 8637-2 specifications.

The performance values should be regarded as approximate.

Results may vary due to differences in measurement method as well as batch variation.

* For additional information see the section under "WARNINGS" in the Instructions for Use.

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
FRESENIUS MEDICAL CARE

Fresenius Medical Care AG
61346 Bad Homburg, Germany
+49 6172 609-0

Haemodialysis 血液透析

en Instructions for Use

FX CorDiax Capillary Dialyser

GENERAL NOTES

Refer to product or carton label for:

	Single use only		Blood pathway sterile. Steam sterilised
	Expiry date		Refer to instructions for use
	Batch		Order number
	Date of manufacture		Storage temperature range +5°C ~ +30°C
	Units		Use only on machines with exact UF-control

Indications: FX CorDiax dialysers are designed for single use in chronic haemodialysis or haemo(dia)filtration respectively.

Contraindications: Special contraindications are unknown. Generally, contraindications for haemodialysis are applicable. Patients with known hypersensitivity to any of the dialyser's material must not be treated with the dialyser.

Side-effects: Certain side effects may occur during dialysis and may result from factors specific to the patient, operating parameters, equipment, priming procedure, dialysis solution, dialyser, anticoagulation, medication etc. Therefore the selection of a dialyser and the selection and monitoring of treatment parameters based on individual patient characteristics, therapy tolerance and clinical requirements as well as compliance with the water and dialysis fluid standards are essential to minimise side effects. The FX CorDiax dialyser is designed for high performance dialysis. In patients not treated with this dialyser before and incident patients starting HD or HDF therapy, treatment intensity shall be gradually increased to permit adequate adaptation. Hypersensitivity or hypersensitivity-like reactions have been observed during dialysis mainly in the first weeks of treatment with the dialyser. Symptomatology can vary and may include: dyspnoea, chest congestion, bronchospasm, respiratory arrest, hypotension, tachycardia, urticaria, erythema, flushing, angioedema, ocular hyperaemia, pruritus, abdominal pain, nausea, convulsions and unconsciousness. Carefully monitor patients who have not previously been treated with the dialyser, or who have shown possible hypersensitivity symptoms during previous treatments, or who have a history of allergy including asthma. If severe hypersensitivity or hypersensitivity-like reactions occur, the dialysis must be discontinued and the blood from the extracorporeal system must not be returned to the patient. Initiate appropriate emergency medical treatment.

The dialyser is steam sterilised and thus contains no sterilisation residues. For patients with a known ethylene oxide hypersensitivity all components of the extracorporeal circuit should be ethylene oxide-free.

Anticoagulation: It is recommended to introduce an anticoagulant to the extracorporeal circuit. Anticoagulation requirements may vary with the patient's condition, application site, dialyser characteristics and treatment modality. Nature, amount and method of application of an anticoagulant must be prescribed by the responsible physician.

Materials: Membrane: Helixone®plus (Polysulfone - PVP blend), Housing: Polypropylene, Potting material: Polyurethane, Sealing-ring: Silicone, Sterility caps: Polypropylene.

Further information may be obtained on request.

WARNINGS

Due to the high hydraulic permeability of the membrane, the dialysers must only be used with dialysis machines which enable precise volumetric fluid control. Refer to dialysis machine manufacturer's technical and safety instructions for the use of high permeability membranes.

During postdilution-HDF and postdilution-HF with manual setting of the filtration flow rate the maximum filtration flow (substitution rate + ultrafiltration rate) must not exceed 25 % of the effective blood flow rate. This is to keep the blood water content within a certain uncritical range and to avoid the risk of haemoconcentration and coagulation in the extracorporeal circuit. The TMP has to be monitored.

Use only if unit package is intact, sealing caps are in place and the dialyser is undamaged. Dialysers must not be used after expiry date (see label). Each dialyser is checked for integrity prior to leaving the factory. If a blood leak should arise, the dialyser must be exchanged.

The dialyser is intended for single-use only. Re-use may be hazardous to both the patient and operator. Cleansing solutions and disinfectants may damage materials employed for the housing, potting and membrane. Safety of use can no longer be guaranteed and the manufacturer assumes no liability. Blood flow rates below the recommended blood flow range can be set. The coagulation should be monitored and the decrease of the clearance values might be steeper than expected.

Blood flow rates above the recommended blood flow range can be set. This might result in a decrease of the sieving coefficient of middle molecular substances (e.g. β_2 -m, myoglobin) and an increase of the sieving coefficient of albumin when performing HDF or HF. The increase of the middle molecule clearance might be lower than expected from raised blood flow and filtration flow rates.

HAEMODIALYSIS

Priming

Affix dialyser in a vertical position. Aseptically connect the arterial and venous blood lines to the lower and upper dialyser ports respectively. Attach the dialysate tubing connectors with the inlet to the upper and the outlet to the lower dialysate port. Fill dialyser according to the instructions provided with the machine (**there is no need to turn the dialyser**). We recommend bag-free priming with ONLINEplus (priming volume: 500 mL). Alternatively the blood compartment can be filled and de-aerated in recirculation mode with isotonic saline solution (e.g. 500 mL bag). Ensure complete de-aeration of the blood compartment.

Patient Connection

Connect arterial line to the patient's circulation. Allow blood to flow in the blood line and dialyser (pump speed approx. 100 mL/min). Connect venous blood line to the patient.

Adjust treatment parameters to desired setting.

The weight loss should be monitored.

Termination

Perform reinfusion of the blood at the end of the treatment according to the instructions provided with the machine. We recommend the bag-free ONLINEplus procedure. Alternatively isotonic saline solution (e.g. 500 mL bag) can be used for reinfusion. The blood should be reinfused completely.

WARRANTY

Products with manufacturing defects will be replaced if the defect is reported with details of the lot number.

The manufacturer will not be liable for any misuse, improper handling, non compliance with instructions for use and cautionary notes and for any damage incurred subsequent to the manufacturer's delivery of the dialyser.

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