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Powder for production of liquid bicarbonate concentrate for bicarbonate haemodialysis

Type of application: For extracorporeal bicarbonate haemodialysis or bicarbonate haemodiafiltration (haemodiafiltration following the FRESenius ONLINE HDF procedure) according to a physician's prescription. Sufficient for a bicarbonate dialysis session of 6h minimum (bibag 650g) and 8h minimum (bibag 900g) at a dialysate flow rate of 500ml/min and a sodium bicarbonate concentration of 32 mmol/L in the ready-to-use haemodialysis solution. For one dialysis session.

Contraindications: Hyperkalaemia (only for haemodialysis concentrates yielding a potassium concentration in the ready-to-use dialysis fluid of 4mmol/L or higher), hypokalaemia (only for haemodialysis concentrates yielding a potassium concentration in the ready-to-use dialysis fluid of 3mmol/L or lower). Potential adverse effects of rapid plasma potassium concentration changes should be considered when prescribing a patient's dialysis fluid potassium concentration. Uncontrollable coagulation anomalies.

Side effects: Hypotension, nausea, vomiting, muscle cramps, dizziness, headache and unconsciousness have been observed during haemodialysis treatment.

Precautions: Only to be used in combination with acid concentrate for bicarbonate haemodialysis in the given dilution according to a physician prescription. Other dilutions can endanger patients. Only to be used with Fresenius Medical Care haemodialysis machines 5008 and/or Fresenius Medical Care haemodialysis machines type 4008 equipped with the adapted connector. For safe use, the dialysis sessions must be carried out by personnel qualified in haemodialysis, conversant with both the current instructions for use and the operating instruction manual of the dialysis machine used.

Note: Acid and basic concentrates for bicarbonate haemodialysis should be diluted immediately prior to use, with water of appropriate quality for haemodialysis. Always check the composition of the ready-to-use haemodialysis solution. Do not use damaged bibag. Discard any residues at the end of dialysis. Failure to do so could result in contamination and cross-infection and could lead to an alteration in the expected performance of the product and pose a risk for the patient for which the manufacturer will assume no responsibility. Powder is non-pyrogenic.

For haemodiafiltration application (FRESenius ONLINE HDF procedure), special attention should be paid to the water quality in order to achieve the requirements of the applicable standards and/or the European Pharmacopoeia for haemodiafiltration solutions.

Storage: Store between +5°C to +30°C. Do not refrigerate. In order to avoid any contamination, keep bibag tightly sealed until it is connected to the haemodialysis machine.

Bag: PA/PE

Presentation: bibag 650g, bibag 900 g

GENERAL INSTRUCTIONS FOR USE FOR bibag DIALYSIS

Observe: More detailed instructions are given in the operating instructions for the haemodialysis machine 5008 and/or Fresenius Medical Care haemodialysis machines type 4008 equipped with the adapted connector.

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Australian Distributor:
Fresenius Medical Care Australia Pty. Ltd
305 Woodpark Road - Smithfield NSW 2164

Reg. No : SIN13822P

DVR-8099
Imported and Distributed by:
Fresenius Medical Care Philippines Inc., 18th Floor Aeon Centre,
corner Alabang-Zapote Road and North Bridgeway, Filinvest
Corporate City, Alabang, Muntinlupa City 1781

KEMENKES RI AKL.20805210485.
Manufactured by Fresenius Medical Care SMAD, France
Imported & distributed by PT. Fresenius Medical Care Indonesia,
Jakarta - Indonesia

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For India Only
Imported By: Fresenius Medical Care India Pvt. Ltd., B-7 & B-8,
Middle Portion, Sipcot Industrial Park, Irungattukottai,
Sriperumbudur-602105, Chennai, India.
Marketed By: Fresenius Medical Care India Pvt. Ltd., S-21, Second
Floor, District Centre, Mayur Palace, Mayur Vihar Phase-1, New
Delhi-110091, India.

In Case of Consumer Complaints, Please Contact: Fresenius
Medical Care India Pvt. Ltd., SAS Tower-A, Ground Floor, The
Medicity, Sector-38, Gurugram-122001, Haryana, India.
Toll Free Number: 18001209500
Email id: customerservice.india@fmc-asia.com
Manufactured By: Fresenius Medical Care SMAD, Z.I. de la
Pontchonnaière, Route de la Chanade / Savigny, 69591 L'Arbresle
Cedex, France.

Shelf life: See Manufacturing Date and Expiry Date on the reverse
of the bag and on the box label.

Broj Rešenja ALIMS: 515-02-02112-16-001
Ime i adresa nosioca upisa:
Fresenius Medical Care Srbija d.o.o.,
Beogradski put bb,
26300 Vršac,
Republika Srbija

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Fresenius Medical Care AG & Co. KGaA